



SaniSure®
Solution Based Innovation

CAP2V8®

BOTTLE SOLUTIONS



From Concept to Prototype to Production

CAP2V8® BOTTLE SOLUTIONS



Table of Contents

PRODUCT CATEGORY	PAGE
STANDARD STERILE MEDIA BOTTLES	4
MEDIA BOTTLE CAPS	10
GL45 GLASS CAPS	12
STANDARD STERILE CENTRIFUGE TUBES	14
CENTRIFUGE CAPS	16
SHAKER FLASKS	18
ACCESSORIES	20

SaniSure® Cap2v8® (say “captivate”), was designed to capture, store and transfer high-purity fluid in a way that is unique to the single-use market. This product was created to solve a problem that you, our customers, had identified as preventing you from successfully making your product. **SaniSure was built around the idea of bringing innovative solutions to our customers — we want to help solve your problems so that you can better achieve your primary goal — to save lives.**

We know that our products are used in critical applications that can directly affect patients, and **everything we do at SaniSure is directed to providing the best possible single-use solutions** to help you create treatments for those patients. That ideal is very important to us and one that we apply to our entire operation.

All Cap2v8 products incorporate the Cap2v8 cap — a one-piece molded tube-and-cap system that provides for aseptic transfer, storage and sampling for **media bottles, glass bottles, centrifuge tubes** and **shaker flasks**. In addition, we provide a wide array of custom bottle solutions to solve even the most demanding problems.

The molded design of the Cap2v8 assembly creates a **closed sterile system that protects product integrity** and provides an aseptic seal for high-purity biologicals, reducing the possibility of contamination and product loss often associated with manual processes and designs. A Cap2v8 bottle assembly reduces your labor and installation time, eliminates cleaning validation and helps you produce a more cost-effective operation.

This catalog highlights some of the standard Cap2v8 products, but these solutions are just the beginning of what Cap2v8 can be. Multiple custom options are available — **from concept to prototype to production**. SaniSure can help you design the perfect solution for your process.

SANISURE — SOLUTION BASED INNOVATION

Sincerely,

Your SaniSure Team

SHAKER FLASKS




STERILE SHAKER FLASKS



Nalgene®

CAP2V8® 45-430 SHAKER FLASKS

CAP2V8 45-430 FOR NALGENE				
CAP2V8	FEED 1 W/ DIP TUBE	FEED 2 W/ DIP TUBE	VENT NO DIP	MATERIAL
45-430 FOR NALGE SHAKER FLASKS	1/8 BY 1/4	NA	1/8 BY 1/4	X,P,A.
45-430 FOR NALGE SHAKER FLASKS	1/4 BY 7/16	NA	1/4 BY 7/16	X,P,A.
45-430 FOR NALGE SHAKER FLASKS	1/4 BY 1/2	NA	1/4 BY 1/2	X,P,A.
45-430 FOR NALGE SHAKER FLASKS	1/4 BY 7/16	NA	1/8 BY 1/4	X,P,A.
45-430 FOR NALGE SHAKER FLASKS	1/8 BY 1/4	1/8 BY 1/4	1/8 BY 1/4	X,P,A.



Corning®

CAP2V8 43 SHAKER FLASKS

CAP2V8 43 FOR CORNING				
CAP2V8	FEED 1 W/ DIP TUBE	FEED 2 W/ DIP TUBE	VENT NO DIP	MATERIAL
43 FOR CORNING SHAKER FLASKS	1/8 BY 1/4	NA	1/8 BY 1/4	X,P,A.
43 FOR CORNING SHAKER FLASKS	1/4 BY 7/16	NA	1/8 BY 1/4	X,P,A.
43 FOR CORNING SHAKER FLASKS	1/4 BY 7/16	NA	1/4 BY 7/16	X,P,A.
43 FOR CORNING SHAKER FLASKS	1/8 BY 1/4	1/8 BY 1/4	1/8 BY 1/4	X,P,A.

*CAP MATERIALS

X = C-Flex® with polypropylene cap
A = AdvantaFlex® with polypropylene cap

P = PharMed® with polypropylene cap
V = Pthalate-free PVC with Hytrel® cap

Contact **SaniSure®** for other solution-based innovations.

SaniSure
Clamp



Single-Use
Tubaset
Assemblies



Mixed4Sure



Transition
Innovative
Solutions



Assessory Kit
for
Glass BioReactor



SaniSure
Gamma Stable
Stir Bar



6914-101918-17