



REVISION HISTORY


| REV | DESCRIPTION | DATE | APPROVED |
|-----|-------------|------|----------|
| | | | |



Replacement Diaphragm for :-

OEM PART = BP2

USED ON 1/4", 3/8", 1/2" (DN4,6,8) Diaphragm Valve

| | | | | |
|--|------------|----------|---|----------|
| | NAME | DATE |  DIAPHRAGM DIRECT | |
| DRAWN | M. MC KEON | 08/21/14 | | |
| CHECKED | | | | |
| <u>Material -</u> <u>EPDM BACKING</u> <u>TFM - DIAPHRAGM</u> | | | TITLE DD03-STFML | |
| Complies with CFR 177.2600 and Conforms to USP Class VI Biological Reactivity Tests Section #87 and #88 | | | SIZE A | REV B |
| | | | FILE NAME: DD03-STFML.dft | |
| SCALE: 1 : 1 | | ECO:0605 | SHEET 1 OF 1 | |

Disclaimer - Diaphragm Direct® does not manufacture, represent, nor sell OEM valves or diaphragms manufactured by any other company. Diaphragm Direct® brand diaphragms are not sponsored by any other manufacturers. Any name, model, part number, brand, image, material or referral is for descriptive purposes only. Diaphragm Direct® brand replacement diaphragms are guaranteed to fit OEM equipment and meet or exceed OEM specifications. Diaphragm Direct® brand replacement diaphragms are covered by the standard Diaphragm Direct® limited product warranty.